



CRITICAL CARE

ACUTE NIV

TRANSPORT



# Flight 60<sup>®</sup>


## Critical and Acute Care Ventilation

The Flight 60 is designed to meet the increasing demand of high quality ventilation while improving patients comfort and independence

### ICU Quality Ventilation Anywhere

- Hospital
  - Acute Care, ER, Transport, Step Down
- EMS and Transport
- Emergency and Preparedness
- Home care

### CCV - Continuity of Care Ventilation



### Full Ventilation Package

- Advanced Ventilation Modes
- Top of the line performance
- User friendly
- Comprehensive Monitoring
- Multipurpose Dual and Single Limb
- High degree of Autonomy - 12 hours internal battery



# FLIGHT 60 Technical Specifications

## Turbine Configuration

| MODES                             |   |
|-----------------------------------|---|
|                                   | ACMV - PCV/VCV/PRVC   |
|                                   | SIMV - PCV/VCV/PRVC   |
|                                   | SPONT (CPAP, BPAP)  |
|                                   | Volume Guarantee  |
|                                   | B-Lev (Bi-Phasic, APRV)                                       |
| CONTROLS                          |   |
| NPPV                              | OFF/LOW/HIGH<br>(leak compensation up to 30 LPM)              |
| VG Mode                           | VtG (Tidal Volume Guarantee)<br>MVG (Minute Volume Guarantee) |
| SIGH                              | ON/OFF  |
| 2min 100% O <sub>2</sub> Function | ON/OFF  |
| Tidal Volume                      | 30 to 2,200 ml  |
| Breath Rate                       | 1 to 99 b/min   |
| Inspiration Time (Ti)             | 0.1 to 3.0 sec  |
| Flow                              | 2 to 100 L/min  |
| Pressure Control                  | 5 to 80 cmH <sub>2</sub> O                                    |
| Pressure Support (PSV)            | 0 to 60 cmH <sub>2</sub> O                                    |
| PEEP/CPAP                         | 0 to 40 cmH <sub>2</sub> O                                    |
| Pressure Trigger                  | -20.0 to -0.1 cmH <sub>2</sub> O                              |
| Flow Trigger                      | 1 to 20 LPM   |
| Rise Profile                      | 5 levels  |
| PSV Ti                            | 0.1 to 3 sec  |
| PSV Flow Termination              | 10% to 90%  |
| Volume Control                    | Time/Flow   |
| Flow Waveform                     | Square/Descending   |
| FiO <sub>2</sub>                  | 21% to 100%   |
| FiO <sub>2</sub> Sensor           | ON, OFF, Calibrate  |
| Manual Breath                     | 0 to 3 sec  |
| Panel Lock                        | ON/OFF  |
| Trends                            | ON/OFF/Clear  |
| Maneuvers                         | Inspiratory Hold, Expiratory Hold                             |
| Hold Length                       | 1 to 6 seconds  |
| VG Mode Controls                  |   |
| Target VtG                        | 30 to 2,200 ml  |
| PSV min                           | 0 to 60 cmH <sub>2</sub> O                                    |
| PSV max                           | 5 to 60 cmH <sub>2</sub> O                                    |
| Trigger Delay                     | ON/OFF  |
| B-Lev Controls                    |   |
| P High                            | 3 to 60 cmH <sub>2</sub> O                                    |
| P Low                             | 0 to 30 cmH <sub>2</sub> O                                    |
| T High (Ti)                       | 1 to 15 sec   |
| T Low                             | 0.5 to 5 sec  |
| Inverse I:E                       | 30:1  |
| ALARMS (variable)                 |   |
| Alarm Prioritization              | 3 Levels – Caution, Medium, High                              |
| Low Minute Volume                 | 0.0 to 50 L/min   |
| Low Pressure                      | OFF, 1 to 98 cmH <sub>2</sub> O                               |
| High Pressure                     | 4 to 99 cmH <sub>2</sub> O                                    |
| High Minute Volume                | 0.1 to 50 L/min   |
| High FiO <sub>2</sub>             | 31% to 99%, OFF   |
| Low FiO <sub>2</sub>              | OFF, 22% to 90%   |
| High Rate                         | OFF, 1 to 99 bpm  |
| Low Rate                          | OFF, 1 to 99 bpm  |
| Low Vte                           | OFF, 10 to 2,200ml  |
| Low Vti                           | OFF, 10 to 2,200ml  |
| Apnea/Back-Up Ventilation         | 10-60 sec   |

| ALARMS (automatic)              |   |
|---------------------------------|---|
|                                 | Check Circuit (Circuit Disconnect), Low/Empty Battery, O <sub>2</sub> Supply Failed, Check O <sub>2</sub> Sensor, Target Volume not reached |
| MONITORED PARAMETERS            |   |
| Waveforms                       | Pressure, Flow, Volume  |
| Loops                           | Pressure/Volume & Flow/Volume   |
| Trends                          | Rate, Peak Inspiratory Pressure & Vte (up to 72 hours)  |
| Airway Pressure LED Gauge       | -10 to 120 cmH <sub>2</sub> O   |
| Peak Inspiratory Pressure       | 0 to 120 cmH <sub>2</sub> O   |
| Base Pressure                   | 0 to 99 cmH <sub>2</sub> O  |
| Mean Pressure                   | 0 to 99 cmH <sub>2</sub> O  |
| Exhaled Tidal Volume            | 0 to 10L  |
| Exhaled Minute Volume           | 0 to 99L  |
| Inhaled Tidal Volume            | 0 to 10L  |
| Inhaled Minute Volume           | 0 to 99L  |
| Actual Breath Rate              | 0 to 99 b/min   |
| Peak Inspiratory Flow           | 1 to 120 L/min  |
| RSBI                            | 0 to 200 1/min*L  |
| Lung Mechanics                  | Static and Dynamic Compliance, Resistance, Plateau Pressure, Auto PEEP  |
| FiO <sub>2</sub>                | 21% to 100%   |
| I:E Ratio                       | 1:99 to 3:1   |
| Battery Level                   | 100% to 0%, Low, Empty  |
| SPECIAL FUNCTIONS               |   |
| Buzzer Level                    | LOW/HIGH  |
| Keypad Buttons                  | Keypad buttons with audible indicator   |
| Power Save                      | ON/OFF/NIGHT  |
| Languages                       | English, French, German, Greek, Hungarian, Italian, Polish, Portuguese, Russian, Spanish, Turkish   |
| Quick Start                     | 5 preset ventilation modes  |
| SIZE AND WEIGHT                 |   |
| Dimensions (WxLxH)              | 29 x 28 x 25 cm / 11.4" x 11.0" x 9.8"  |
| Weight                          | 6.3 Kg / 13.9 lbs<br>6.9 Kg / 15.2 lbs (with integrated mixer)  |
| OXYGEN                          |   |
| O <sub>2</sub> Mixer (optional) | Internal integral, electronically controlled  |
| High Pressure                   | 35 to 90 psi  |
| Low Flow Port                   | 0 to 15 L/min   |
| Low Flow Blending Bag           | 0 to 15 L/min   |
| POWER SUPPLY                    |   |
| AC Power Inlet                  | 100 to 240 VAC, 50-60Hz   |
| DC Power Inlet                  | 12 to 15 VDC  |
| Internal Batteries              | Hot Swappable 12 hours operation  |
| Charging Time                   | Up to 3 hours   |
| COMMUNICATIONS                  |   |
| USB x2                          | Download Logs, SW Upgrade   |
| RS232 x2                        | Remote Alarm and Monitoring   |
| ENVIRONMENTAL                   |   |
| Operation Temperature           | -18°C to 50°C / -0.4°F to 122°F   |
| Storage Temperature             | -20°C to 71°C / -4.0°F to 160°F   |
| Relative Humidity               | 15% to 95% at 31°C / 88°F   |
| Operation Altitude              | 110 kPa to 70 kPa   |
| Water/Dust Resistance           | IP34 (splash proof)   |
| STANDARDS                       |   |
|                                 | IEC 60601-1, IEC 60601-1-2, IEC 60601-2-12, ASTM 1246F, ISO 10651-2/3, RTCA DO-160 F  |

\* NOTE: Some features of the Flight 60 may not be available or cleared in all countries. For detailed information, contact your local Flight Medical representative.